

Comparison of a *Histoplasma/Blastomyces* Lateral Flow Assay to Enzyme Immunoassay Using Dog and Cat Sera

Janelle Renschler, Logan Cheney, Amanda Albers, Michelle Durkin and L. Joseph Wheat

MiraVista Diagnostics, Indianapolis, Indiana, USA

INTRODUCTION

- Quantitative enzyme immunoassays (EIAs) for *Histoplasma* and *Blastomyces* antigen are highly sensitive for the detection of fungal antigen to aid in the diagnosis of histoplasmosis and blastomycosis in dogs and cats.¹⁻³ These assays show nearly complete cross reactivity; therefore, both infections may be detected by a single antigen test.
- A serum lateral flow assay (LFA) was developed for human histoplasmosis, and this might be used in dogs and cats.
- Use of the LFA as a point-of-care assay in veterinary clinics would allow for more rapid diagnosis and earlier treatment of systemic mycoses
- The objective of this study was to compare agreement between the *Histoplasma* antigen LFA and the quantitative EIAs for *Histoplasma* and *Blastomyces*, using canine and feline serum.

METHODS

- Serum samples previously tested in either EIA, each tested in duplicate by LFA
 - 11 EIA-negative (22 LFA devices)
 - 21 EIA-positive (42 devices)
 - Ranged from 0.2 ng/mL to 19 ng/mL
 - Species distribution: 12 feline, 20 canine
- Serum was treated with EDTA and extracted supernatant was loaded onto the LFA device
- LFA devices were interpreted:
 - Visually by three qualified human evaluators (Fig. 1)
 - Two automated readers (Cube-Reader, opTricon GmbH, Germany) (Fig. 2).
- A **consensus RANN score** (ranging from 0-10 based on intensity of the test line compared to scorecard standards) was obtained by using two evaluator agreement. The samples were randomized and evaluators were blinded to EIA results.
- Cube reader T_{mean} also recorded and reported as positive or negative based on cutoff calculation

RESULTS

- Results comparison summarized in Table 1
- Kappa agreement with EIA:
 - Visual assessment **0.858**
 - Automated readers **0.861** and **0.760**
- EIA-negative samples with positive LFA results included:
 - Three devices with very low consensus RANN scores (1-2)
 - Two were also positive by both cube readers
 - One was positive by only one cube reader
 - All were canine samples
- EIA-positive samples with negative LFA results included:
 - One device with very low positive EIA result (0.5 ng/mL) and negative results by all LFA evaluations
 - Three additional devices with very low positive EIA results, negative by one cube reader
 - All were feline samples



Fig. 1. LFA devices with RANN scores ranging from negative to high positive



Fig. 2. opTricon Cube-Reader placed over the LFA device

Table 1. Comparison of *Histoplasma* or *Blastomyces* EIA results with LFA by visual assessment or opTricon Cube-Readers

EIA Result	LFA Visual Result		LFA Cube 1 Result		LFA Cube 2 Result	
	-	+	-	+	-	+
-	19	3	20	2	19	3
+	1	41	2	40	4	38

CONCLUSIONS

- Good agreement between LFA and EIA
- Uncommon false positive LFA results from some canine samples; very faint bands
 - also detected by automated readers (not evaluator error)
- Rare false negative LFA results from samples with low antigen concentration
 - Visual assessment outperformed automated readers
- LFA may be a useful diagnostic tool for rapid serum diagnosis of histoplasmosis or blastomycosis in dogs and cats**

MiraVista  **DIAGNOSTICS**
Rapid Fungal Testing. Accurate Results.

References:

- Spector D, Legendre AM, Wheat J, et al. Antigen and antibody testing for the diagnosis of blastomycosis in dogs. *J Vet Intern Med* 2008;22:839-43.
- Cook AK, Cunningham LY, Cowell AK, et al. Clinical evaluation of urine *Histoplasma capsulatum* antigen measurement in cats with suspected disseminated histoplasmosis. *J Feline Med Surg* 2012;14(8):512-5.
- Cunningham L, Cook A, Hanzlicek A, et al. Sensitivity and specificity of *Histoplasma* antigen detection by enzyme immunoassay. *J Am Anim Hosp Assoc*. 2015;51(5):306-10.